




































| | | |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| | | |
| SC-01 | Upper cleaning chassis,without spray arm, place the module cleaning basket |  |
| SC-02 | Upper cleaning chassis,with spray arm, place the module cleaning basket |  |
| XC-03 | Lower cleaning chassis,without spray arm, place the module cleaning basket |  |
| PLYT-08 | Spray bottle module basket 8 positions Bottle size: 250~2000ml Used to clean blue cap bottle, erlenmeyer flask, volumetric flask, triangular flask, flask,etc.Equipped with three-petal flower basket bottle holder of stainless steel and cross bottle holder. Injection nozzle:Φ8×H165mm*8 Designed according to blue cap bottle of 1000 ml, each module can clean 21 blue cap bottle. |  |
| PLYT-15 | Spray bottle module basket 15 positions Bottle size: 50~250ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask,equipped with three-petal flower basket bottle holder of stainless steel and cross bottle holder. Injection nozzle:Φ4.6×H165mm*21 Designed according to volumetric flask of 250 ml, each module can clean 21 volumetric flask. |  |
| PLYT-24 | Spray bottle module basket 24 positions Bottle size: 50~250ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask,equipped with three-petal flower basket bottle holder of stainless steel and cross bottle holder. Injection nozzle:Φ4.6×H165mm*32 Designed according to volumetric flask of 100 ml, each module can clean 32 volumetric flask. |  |

| | | |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <p>PLYT-35</p> | <p>Spray bottle module basket 35 positions Bottle size: 10~100ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with three-petal flower basket bottle holder of stainless steel and cross bottle holder. Injection nozzle: $\Phi 3.2 \times H110\text{mm} \times 50$ Designed according to volumetric flask of 50 ml, each module can clean 50 volumetric flask.</p> |  |
| <p>PLYT-48</p> | <p>Spray bottle module basket 48 positions Bottle size: 10~50ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with three-petal flower basket bottle holder of stainless steel and cross bottle holder. Injection nozzle: $\Phi 3.2 \times H110\text{mm} \times 66$ Designed according to volumetric flask of 25 ml, each module can clean 66 volumetric flask.</p> |  |
| <p>PLYT-25</p> | <p>Sample bottle basket 25 positions Bottle size: 500ml, equipped with cross bottle holder Injection nozzle: $\Phi 4.6 \times H165\text{mm} \times 25$ Designed according sampling bottles of 500ml, only one basket can be placed on each layer, and each basket can clean 25 sampling bottles of 500ml.</p> |  |
| <p>PLYT-20</p> | <p>Sample bottle basket 20 positions Bottle size: 1000ml, equipped with cross bottle holder Injection nozzle: $\Phi 4.6 \times H165\text{mm} \times 20$ Designed according to sampling bottles of 1000ml, only one basket can be placed on each layer, and each basket can clean 20 sampling bottles of 1000ml</p> |  |
| <p>PLWT-15</p> | <p>Spray bottle module basket 15 positions Bottle size: 50~250ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with cross bottle holder. Injection nozzle: $\Phi 4.6 \times H165\text{mm} \times 18$ Designed according to volumetric flask of 250 ml, each module can clean 18 volumetric flask.</p> |  |
| <p>PLWT-24</p> | <p>Spray bottle module basket 24 positions Bottle size: 50~250ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with cross bottle holder. Injection nozzle: $\Phi 4.6 \times H165\text{mm} \times 32$ Designed according to volumetric flask of 100 ml, each module can clean 32 volumetric flask.</p> |  |

| | | |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <p>PLWT-48</p> | <p>Spray bottle module basket 48 positions Bottle size: 10~50ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with cross bottle holder. Injection nozzle: $\Phi 3.2 \times H110\text{mm} \times 66$ Designed according to volumetric flask of 25 ml, each module can clean 66 volumetric flask.</p> |  |
| <p>PLWT-35</p> | <p>Spray bottle module basket 35 positions Bottle size: 10~100ml Used to clean erlenmeyer flask, volumetric flask, triangular flask, flask, equipped with cross bottle holder. Injection nozzle: $\Phi 3.2 \times H110\text{mm} \times 50$ Designed according to volumetric flask of 50 ml, each module can clean 50 volumetric flask.</p> |  |
| <p>PLWT-66-1</p> | <p>Spray bottle module basket 66 positions Bottle size: 10~20ml Used to clean test tube, equipped with cross bottle holder.</p> |  |
| <p>PLWT-66-2</p> | <p>Spray bottle module basket 66 positions Bottle size: 10~50ml Used to clean pipette</p> |  |
| <p>PL-150</p> | <p>Spray bottle module basket 150 positions Bottle size: 2ml Used to clean the chromatographic injection bottle Injection nozzle: $\Phi 2 \times H45\text{mm} \times 150$</p> |  |
| <p>PL-102-1</p> | <p>Spray bottle module basket 102 positions Used to clean the chromatographic injection bottle Injection nozzle: $\Phi 2 \times H45\text{mm} \times 102$</p> |  |
| <p>PL-102-2</p> | <p>Spray bottle module basket 102 positions Bottle type: test tube For cleaning test tubes, centrifuge tubes</p> |  |

| | | |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <p>PL-102-3</p> | <p>Spray bottle module basket 102 positions Bottle type: pipettes, graduated pipettes For cleaning pipettes, graduated pipettes.</p> |  |
| <p>MX-102-1</p> | <p>Hybrid module 102 positions Bottle size: 0.1~50ml Used to clean 78 injection bottles and 24 pipettes Injection nozzle: $\Phi 2 \times H 45 \text{mm} \times 102$</p> |  |
| <p>MX-102-2</p> | <p>Hybrid module 102 positions Bottle size: 0.1~50ml Used to clean 78 injection bottles and 24 pipettes</p> |  |
| <p>MX-102-3</p> | <p>Hybrid module 102 positions Bottle size: 0.1~50ml Used to clean 54 injection bottles and 48 pipettes</p> |  |
| <p>MX-102-4</p> | <p>Hybrid module 102 positions Bottle size: 0.1~50ml Used to clean 54 injection bottles and 48 test tube</p> |  |

| | | |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <p>MX-102-5</p> | <p>Hybrid module 102positions Bottle size: 0.1~50ml Used to clean 78 test tube and 24 pipette</p> |  |
| <p>MX-102-6</p> | <p>Hybrid module 102positions Bottle size: 0.1~50ml Used to clean 54 test tube and 48 pipette</p> |  |
| <p>MX-48-1</p> | <p>Hybrid module 48 positions Bottle size: 0.1~100ml Used to clean 24 pipettes and 24 volumetric flasks of 100ml Injection nozzle:Φ4×H165mm*24</p> |  |
| <p>MX-48-2</p> | <p>Hybrid module 48 positions Bottle size: 0.1~100ml Used to clean 24 pipettes and 24 volumetric flasks of 100ml Injection nozzle:Φ4×H165mm*24</p> |  |
| <p>MX-74-1</p> | <p>Hybrid module 74 positions Bottle size: 2~100ml Used to clean 50 injection bottles and 24 volumetric flasks of 100ml Injection nozzle:Φ4×H165mm*24 Injection nozzle:Φ2×H45mm*50</p> |  |
| <p>MX-74-1</p> | <p>Hybrid module 74 positions Bottle size: 2~100ml Used to clean 50 injection bottles and 24 volumetric flasks of 100ml Injection nozzle:Φ4×H165mm*24 Injection nozzle:Φ2×H45mm*50</p> |  |

| | | |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <p>K-02</p> | <p>1/2 Cleaning frame module large quantity Bottle size: 10~2000ml Used to clean jars, glass rods, bottle caps Two 1/2 modules form a layer of cleaning basket</p> |  |
| <p>K-03</p> | <p>1/4 Fan-shaped test tube cleaning frame module with spring Bottle size: 10~20ml Used to clean test tubes</p> |  |
| <p>K-04</p> | <p>1/4 Fan-shaped test tube cleaning frame module large quantity Bottle size: 10~20ml Used to clean test tubes, centrifuge tubes, etc.</p> |  |
| <p>K-01</p> | <p>Cleaning frame module large quantity Bottle size: 10~2000ml Used to clean jars, glass rods, bottle caps Only one basket can be placed on each layer</p> |  |
| <p>K-30</p> | <p>1/2 Cleaning frame module 30 positions Bottle size: petri dish with a diameter of 60~150mm Each module can wash 30 petri dishes of 90mm Two 1/2 modules form a layer of cleaning basket</p> |  |